# Sonopuls 190

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo



Copyright:

**Enraf  
Nonius**

Enraf-Nonius B.V.  
P.O. Box 810  
2600 AV DELFT  
The Netherlands  
Tel.: +31(0)15 - 26 98 400  
Fax: +31(0)15 - 25 62 686

Order number : 1416.750-41  
1<sup>st</sup> print, December 2000

# Sonopuls 190

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

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## TABLE OF CONTENTS

<b>1. Introduction</b>	<b>17</b>
1.1. General	17
1.2. Therapeutic possibilities	17
1.3. Treatment head	17
1.4. Efficiency and safety	17
1.5. Programs	17
1.6. Finally	17
<b>2. Preliminary notes</b>	<b>17</b>
<b>3. Installation and safety</b>	<b>17</b>
3.1. Installation	17
3.2. Connection	17
3.3. Technical maintenance	18
3.4. Electromagnetic interference	18
3.5. Product Liability	18
<b>4. The therapy possibilities in a nutshell</b>	<b>18</b>
4.1. Ultrasound therapy	18
<b>5. Indications and Contra-indications</b>	<b>19</b>
5.1. Indications	19
5.2. Absolute specific contra-indications	19
5.3. Relative specific contra-indications	19
<b>6. Controls</b>	<b>19</b>
6.1. APPARATUS	19
<b>7. Installation</b>	<b>20</b>
7.1. Connection of the mains adapter	20
7.2. Switching on and self test	20
<b>8. Disconnection of mains supply</b>	<b>20</b>
<b>9. Operation</b>	<b>21</b>
9.1. Set up the Sonopuls by hand	21
9.2. Pre-programmed protocols selection	22
9.3. Overview treatment suggestions	22
9.4. Programming	23
9.5. Contact control; sound signal on/of	23
9.6. Actions in relation to treatments	24
<b>10. Remarks on operation</b>	<b>25</b>
10.1. The treatment head: a precision instrument	25
10.2. The contact medium	25
10.3. Contact control at incomplete contact	25
<b>11. Maintenance by the user</b>	<b>25</b>
11.1. Apparatus	25
11.2. Treatment head	25
11.3. Environment Information	25
<b>12. Fault conditions</b>	<b>25</b>
12.1. Display's fail to light up	25
12.2. Error code ErX* on the duty cycle display	25
12.3. The contact control fails to operate	25
<b>13. Specifications</b>	<b>26</b>
13.1. Ultrasound	26
13.2. Treatment head	26
<b>14. Technical data</b>	<b>26</b>
14.1. Mains adapter	26
14.2. Sonopuls	26
14.3. Environment conditions for Transport and Storage	27
14.4. Environment conditions normal use	27
<b>15. Ordering data</b>	<b>27</b>

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## 1. INTRODUCTION

### 1.1. General

The Sonopuls 190 is a unit for continuous and pulsed ultrasound therapy at a set frequency of 1 MHz. All functions of the unit are controlled and monitored by a built-in microprocessor, ensuring a high degree of reliability and safety.

### 1.2. Therapeutic possibilities

By using continuous and pulsed ultrasound it is possible to achieve muscle relaxation and pain reduction and also to improve the blood perfusion and recovery ability of tissues.

The Sonopuls 190 cannot be connected to an electrotherapy unit. It is therefore not possible to apply combinationtherapy with the Sonopuls190.

### 1.3. Treatment head

The large treatment head is permanent connected to the unit. The head has excellent beam characteristics, fully meeting the requirements of the existing standards. The treatment head is suitable for subaquatic treatments. The connection of a second, small treatment head, is not possible.

### 1.4. Efficiency and safety

During treatment the Sonopuls 190 constantly checks the contact control. As soon as the acoustic contact is poor the therapist is warned by means of an audio signal (if turned on, see chapter "operation"), the intensity will be reduced automatically and the treatment timer will stop.

The contact control in combination with the excellent beam characteristics, maintains the quality of the treatment. Maximum safety is inherent in the design of the unit. Both the electrical safety and the ultrasound safety of the Sonopuls 190 comply fully with the IEC safety standards.

### 1.5. Programs

In the Sonopuls 190 9 easy adjustable treatment suggestions are available which serve as a guideline for treatments. A free programmable start up program ensures that, after switching on, the unit automatically sets itself to the most frequently used therapy type. For setting of the start up program and an overview of the pre-programmed programs we refer to chapter 9 "Operation" of this manual.

### 1.6. Finally

You have made a wise choice in selecting the Sonopuls 190. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorised distributor.

## 2. PRELIMINARY NOTES

It is important that you read these operating instructions carefully before using the Sonopuls 190. Please make sure that these instructions are available to all personnel who operate the equipment.

For the various applications of this apparatus, please refer to our therapy book "Ultrasound Therapy" reference number 1482.762.

Pay attention to the following before using the Sonopuls 190:

1. Keep yourself informed of the contra-indications (see page 19).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to short-wave or microwave equipment.
3. The apparatus may not be used in so-called "wet-rooms" (hydrotherapy rooms).

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

## 3. INSTALLATION AND SAFETY

### 3.1. Installation

- Do not install the unit in a location near to a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- Cordless telephones should not be used in close proximity of the unit.
- This unit should not be used in so-called "wet-rooms" (hydrotherapy rooms). The unit has to be installed in such way that liquid cannot enter.
- Always use the by Enraf-Nonius supplied mains adapter (type ENB-1530).

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket (if connected) and have it checked by an authorised expert.

### 3.2. Connection

- Mains supply connections must comply with the national requirements regarding medical rooms.
- Before connection of this apparatus to the supply, check that the mains voltage and frequency stated on the type plate of the mains adapter corresponds to that of the mains supply.
- The mains adapter is a part of the supply circuit on which the safety of the equipment is based. As a consequence, the certificates (CB and GOST) of the Sonopuls 190 are only valid in combination with this particular type of mains adapter (ENB-1530).

**It is not permitted to connect another mains adapter than the type prescribed by Enraf-Nonius for the Sonopuls 190.**

### 3.3. Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by another agency, authorised by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

Maintenance and all repairs should only be carried out by an authorised agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorised persons.

Opening of the equipment by unauthorised agencies is not allowed and will terminate any claim to warranty.

### 3.4. Electromagnetic interference

Operation in close proximity (i.e. less than 2 metres) to working short-wave or microwave therapy units may produce instability in the Sonopuls 190 output. To prevent electromagnetic interference, we strongly advise using separate mains groups (phases) for the Sonopuls 190 and the short-wave or microwave equipment, and a distance of at least 2 metres is kept between the Sonopuls 190 and the short-wave or microwave equipment. Ensure that the mains cable of the short-wave/microwave unit does not come near the Sonopuls 190 or the patient.

If problems with electromagnetic interference persist, please contact your authorised distributor.

### 3.5. Product Liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

#### U.S.A. only:

##### Safety Precautions

This unit operates with high voltages. Servicing of the Sonopuls 190 should be referred to qualified service technicians or returned directly to the distributor. To assure continued compliance with FDA 21 C.F.R. 1050.10 standards, the Sonopuls 190 should be calibrated and safety tested at least once each year. The producer recommends that all service work be performed by the distributor. However, a service manual can be obtained by writing to the factory for those who are qualified to service the device and are fully aware of their responsibilities.

##### Caution:

Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.

## 4. THE THERAPY POSSIBILITIES IN A NUTSHELL

### 4.1. Ultrasound therapy

#### Frequency

With the Sonopuls 190 ultrasound therapy with a ultrasound frequency of 1 MHz can be applied.

#### Ultrasoundmode

Both continuous and pulsed ultrasound can be applied. The advantage of the pulsing sound beam is that the thermal sensation is suppressed. In this way it is possible to set a higher intensity than in the case of continuous ultrasound.

#### Pulse frequency / duty cycle

The pulse frequency is set at 100 Hz. The duty cycle (= ratio pulse duration - pulse interval) can be altered in 5, 10, 20, 50, 80 and 100 % (100% = continuous).

## 5. INDICATIONS AND CONTRA-INDICATIONS

Please refer to the optional therapy book.

### 5.1. Indications

- conditions affecting joints, bone, muscle tissue;
- rheumatoid arthritis (not acute);
- affections of the peripheral nerves;
- affections of the blood circulation;
- affections of the internal organs;
- affections of the skin, scar tissue;
- Dupuytren's contracture;
- open wounds, decubitus ulcer, post-traumatic conditions.

### 5.2. Absolute specific contra-indications

- eyes;
- heart;
- pregnancy;
- epiphysial discs;
- brain tissue;
- testicles.

### 5.3. Relative specific contra-indications

- status post laminectomy;
- loss of sensation;
- endoprostheses;
- tumours;
- post-traumatic sequelae;
- thrombophlebitis and varices;
- septic inflammation;
- diabetes mellitus.

## 6. CONTROLS

### 6.1. APPARATUS

(see fold-out sheet at front)

#### [1] Connector for mains adapter

For the connection of the mains adapter.

#### [2] On/Off button

With this button the Sonopuls 190 is turned on or off.

#### [3] Treatment time display

When there is adequate ultrasound contact, the timer will count down to zero. If the ultrasound contact becomes inadequate, the timer stops automatically (a flashing dot in the display indicates that the clock is running).

#### [4] Treatment time buttons

These buttons are used to set the treatment time. Treatment times below 10 minutes can be set in intervals of half a minute by using the

 button.

#### [5] Intensity display

Depending on the choice of display function, this display shows intensity (W/cm<sup>2</sup>) or the power (W). The setting flashes and/or a sound signal is audible (if turned on) as long as there is no adequate acoustic ultrasound contact (contact control).

#### [6] Display function button (W/cm<sup>2</sup> or W)

Switches the display function between intensity (W/cm<sup>2</sup>) and power (W).

#### [7] Intensity buttons

These  /  buttons are used to set the intensity. The maximum intensity can be set from zero by one touch on the  button.

#### [8] Program button

With this button the selection of the program modus is selected. To select one of the pre-programmed suggestions press one of the buttons [7] within 2½ seconds. The indication P1 - P9 is visible in the display [5]. After approx. 2 seconds the treatment parameters are displayed and treatment can be started.

#### [9] Display duty cycle ultrasound and program number

Display of the set duty cycle (= ratio pulse duration - pulse interval) and program number.

#### [10] Duty cycle button

This button is used to set the duty cycle (5, 10, 20, 50, 80 and 100 %).

---

## 7. INSTALLATION

---

### 7.1. Connection of the mains adapter

- Connect the supplied mains adapter to the connector [1].
- Connect the mains adapter to a wall socket. The flashing dot in the upper display [9] indicates that the unit is connected to the mains supply and stand-by.

### 7.2. Switching on and self test

- Switch on the unit using the On/Off button [2].
- Immediately after switching on, the unit carries out a self test.

#### Check that:

- All LED's light up briefly.
- At the end of the test a beep can be heard.

Contact your supplier if this is not the case.

---

## 8. DISCONNECTION OF MAINS SUPPLY

---

- Switch the Sonopuls 190 off with the on/off button [2].
- Unplug the mains adapter from the wall socket.
- Disconnect the mains adapter from the Sonopuls 190.

**9.1. Set up the Sonopuls by hand****Switching on**

- Connect the supplied mains adapter
- Connect the mains adapter to a wall socket
- Switch on the unit using the On/Off button

Immediately after switching on, the unit carries out a self test.

**Choice of therapy mode**

- Select the desired duty cycle  
continuous ultrasound is 100 %, pulsed ultrasound is 5, 10, 20, 50 and 80 %.

The selection is indicated on the display.

**Treatment time**

- Select the treatment time.

The set time is indicated on the display.

**Intensity**

- Select display in watt or W/cm<sup>2</sup>
- Set the intensity /power (the maximum intensity which can be set for a duty cycle of 80 and 100 % is 2 W/cm<sup>2</sup>)

The display will flash until adequate contact is made.

**Starting treatment**

- Position the treatment head (do not forget to apply gel)
- When adequate contact is made, the timer will start (flashing dot)

**Ending treatment**

Treatment can be stopped by taking off the treatment head from the patient and setting the timer to zero. Treatment stops automatically when the treatment time has elapsed.

## 9.2. Pre-programmed protocols selection

Enraf-Nonius stored parameters on the P1 until P9. The parameters can be changed see paragraph 9.4.

### Switch on


- Connect the standard main adapter to connector [1]
- Connect the main adapter plug to the mains socket
- Switch the unit on with the button [2]

During the start-up period runs a self-test.

### Program selection

- Push the button [8]
- Select a program P1 - P9 with one of the buttons [7] within 2½ seconds. The number is indicated on display [5]
- After finishing selections with the buttons [7] the parameters will show after approx. 2 seconds

### Starting treatment

- Position the treatment head on the spot that has to be treated (do not forget to apply gel)
- When adequate contact is made, the timer will start (flashing dot)

### Ending treatment

Treatment can be stopped by taking off the treatment head from the patient and setting the timer to zero. Treatment stops automatically when the treatment time has elapsed.

## 9.3. Overview treatment suggestions

Program	Indication	Duty Cycle	Treat. time [min.]	Intensity [W/cm²]
P1	Arthritis	100%	10	1.0
P2	Bursitis / Capsulitis	10%	8	1.0
P3	Epicondylitis	20%	4	0.5
P4	Contusions / Distorsions	20%	6	1.0
P5	Ligament sprain	50%	6	1.5
P6	Scar tissue	20%	8	0.8
P7	Muscle strain	80%	5	0.4
P8	Spondylosis	20%	10	1.0
P9	Tendinitis	50%	5	0.4

### Attention:

- During treatment the patient may not feel unpleasant sensations amounting to pain. A mild sensation of excitation is permissible.
- If, as result of treatment, headache, vertigo, fatigue and/or other (autonomic nervous) reactions develop, subsequent treatment should be given at a lower intensity.
- With continuous and pulsed ultrasound at high intensity a sensation of heat may be felt. Only a mild sensation of warmth is acceptable.

## 9.4. Programming

The program P0 can be programmed. This position is ideal to store the most frequent used parameters. When the unit start-up these parameters of program P0 will be displayed. The program numbers of the therapy suggestions P1 – P9 can be overwritten the same way.


### Switch on

- Connect the standard main adapter to connector [1]
- Connect the main adapter plug to the mains socket
- Switch the unit on with the button [2]

### Program selection

- Push the button [8]
- Start Selecting the program you wish to change P0 - P9 with one of the buttons [7] within 2½ seconds
- After finishing selections with the buttons [7] the parameters will show after approx. 2 seconds

### Treatment parameters: Adjusting or changing

- Adjust the parameters as decrypted under 9.1.  
Set up the Sonopuls by hand” (follow the steps until adjusting the intensity)

### Storing

- Push at the program button [8] during minimal 3 seconds
- The parameter will be stored under the selected program number. The storing will be confirmed by a double beep


## 9.5. Contact control; sound signal on/of

Inadequate acoustic contact can be recognised by a sound signal; a short beep every 2½ seconds. As desired this sound signal can be turned on or off.

### Sound signal on/off

- Switch the unit on with the On/Off button
- Switch the unit ON while holding down the “duty cycle button” [10] at the same time
- Release the button [10] as soon as in the upper display a “half 8” appears

In this way the sound signal is alternately turned on or off \*(default setting is “on”).

## 9.6. Actions in relation to treatments

### Before treatment

- Check the patient for absolute and relative contraindications
- Test the thermal sensitivity of the treatment area
- The skin of the area in question is cleaned (removal of grease) with soap or 70 % alcohol to permit optimal transmission of the ultrasound
- If the skin is hairy shaving is best

### During treatment

- The treatment head is kept in slow continuous motion, also for the semi-static method
- The patient is regularly asked to report any sensations felt. If necessary the treatment is modified; the intensity can be reduced or a switch is made from continuous to pulsed ultrasound
- If there are indications of poor transfer of the ultrasound energy the contact medium can be reviewed if necessary or it can be moved with the treatment head

### After treatment

- The patient's skin and the treatment head are cleaned with a towel or tissue. The treatment head is also cleaned with 70 % alcohol
- The expected effects are checked (e.g. pain, circulation and mobility)
- The patient is asked to comment subsequently on any reactions that may occur

## 10. REMARKS ON OPERATION

### 10.1. The treatment head: a precision instrument

A treatment head is a precision instrument. Great care is taken in the development and production in order to obtain the best possible beam characteristics. Rough treatment (jarring or dropping) can adversely affect these characteristics, and must therefore be avoided.

### 10.2. The contact medium

In order to ensure efficient transfer of energy, a contact medium is required between the treatment head and the body. Air causes virtually total reflection of the ultrasound energy. The best medium for the transfer of ultrasound energy is a gel.

- For preference, use Enraf-Nonius Contact-Gel®, as this allows the excellent characteristics of the treatment heads to be used to their full advantage.
- The gel should be applied to the part of the body to be treated and then spread out with the treatment head.
- Never apply the gel to the treatment head. The treatment head will register this as contact and may emit ultrasound energy, which could damage the treatment head.

If the body surface is very irregular, making it difficult to obtain good contact between the treatment head and the body, or if direct contact must be avoided (e.g. due to pain), the affected area may be treated under water (subaqual method). The water should be degassed (by previous boiling) in order to prevent air bubbles arising on the treatment head and the body.

### 10.3. Contact control at incomplete contact

When the treatment head only contacts partially with the body, the treatment is interrupted. When not enough surface of the treatment head makes contact with the patient the treatment is interrupted (= contact control). The intensity display starts flashing (and/or a sound signal is audible if turned on) and the treatment time is interrupted. The contact control does not function at a set intensity of 0,1 and 0,2 W/cm<sup>2</sup>.

## 11. MAINTENANCE BY THE USER

### 11.1. Apparatus

Before performing any user maintenance, switch off the apparatus and disconnect the plug from the mains supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

### 11.2. Treatment head

The treatment head and cables should be regularly inspected for damage, e.g. hairline cracks, which could allow penetration by liquids. Clean the contact surface immediately after each treatment. Make sure that no ultrasound gel remains on the treatment head. We further recommend cleaning the head and cable daily, using lukewarm water. The treatment head can be disinfected using a cloth moistened with 70% alcohol.

### 11.3. Environment Information

Your Sonopuls 190 contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart the Sonopuls 190 and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

**Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.**

## 12. FAULT CONDITIONS

### 12.1. Display's fail to light up

Check if the mains adapter is connected to the unit and mains supply.

### 12.2. Error code ErX\* on the duty cycle display

The unit has discovered a fault using its "self-test". Switch the unit back on. If the code re-appears, contact your supplier. The unit is probably defective.

\* ErX: number 1/9 is mentioned on position X

### 12.3. The contact control fails to operate

- Check the contact control by placing the treatment head in a bowl of water. It may be that the contact control is working, but that the contact gel is inadequate.
- An intensity of 0.1 or 0.2 W/cm<sup>2</sup> is used.

**13.1. Ultrasound**

Frequencies	: 1 MHz
Display	: intensity in W/cm <sup>2</sup> and power output in W (SATP*)
Contact control threshold	: 65 %
Timer	: 0 - 30 minutes, coupled to contact control

**Ultrasound, continuous**

Pulse frequency / duty cycle	: 100 Hz / 100 %
Intensity	: 0 - 2 W/cm <sup>2</sup>

**Ultrasound, pulsed**

Pulse frequency / duty cycle	: 100 Hz / 5, 10, 20, 50, 80 %
Intensity	: 0 - 3 W/cm <sup>2</sup> , duty cycle 5, 10, 20, 50 %
	: 0 - 2 W/cm <sup>2</sup> , duty cycle 80 %

**13.2. Treatment head**

1 MHz, large	: Geometric surface 6,2 cm <sup>2</sup> , ERA** 5,0 cm <sup>2</sup> , BNR*** max. 6 type collimating, side panel radiation max. 10 mW/cm <sup>2</sup>
--------------	---

\* SATP = Spatial Average Temporal Peak (average pulse power)

\*\* ERA = Effective Radiating Area, this is the effective radiating area of the treatment head

\*\*\*BNR = Beam Nonuniformity Ratio, indicates the ratio between the peaks and the average value of the intensity in the sound beam. A low BNR excludes high, unwanted energy concentrations.

## 14. TECHNICAL DATA

**14.1. Mains adapter**

Model number	: ENB-1530
Mains voltage	: 100 - 240 Volt
Frequency	: 50/60 Hz
Max, Power consumption	: 40 VA

**14.2. Sonopuls**

Patient leakage current	: typically 1 $\mu$ A (requirement IEC $\leq$ 100 $\mu$ A)
idem single fault condition	: typically 2 $\mu$ A (requirement IEC $\leq$ 500 $\mu$ A)
Safety class	: II* type BF**, according to IEC 60601-1
Dimensions	: 220 x 200 x 195 mm
Weight	: 1.7 kg

Certificates	: CB / GOST
CE marking	:  ***
Class	: IIb ***
Safety standards	: IEC 60601-1 and IEC 60601-2-5.

### 14.3. Environment conditions for Transport and Storage

Environment temperature	: -10° till +50° C
Relative humidity	: 10 till 100 %
Atmospheric pressure	: 500 till 1060 hPa

### 14.4. Environment conditions normal use

Environment temperature	: 10° till 40° C
Relative humidity	: 10 till 90 % (non condensing)
Atmospheric pressure	: 500 till 1060 hPa

\* II : Safety class II (double insulated)

\*\* BF : The equipment has a floating patient circuit.

Type BF

\*\*\*  : According to European requirement 93/42 EEC

*Technical modifications reserved.*

---

## 15. ORDERING DATA

---

For the ordering data of the Sonopuls 190, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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